

	Platform	Programs	Music Program	Autophone	Fitting Channels	Bandwidth (Hz)	Identities/Rationales	Automatic Adaptation Manager	Floating Linear Gain	Directionality Bands	Directionality Modes	Front Focus	Noise Management	Feedback Management	Spatial Sound Compression	Spatial Noise Management	My Voice	Life Learning	Binaural Processing	Binaural Synchronisation	Binaural Coordination	Battery Size	Telecoil	Program Button	Volume Control	ConnectLine Enabled	nEARcom	POWER
Oticon	Agil Pro	4	✓	✓	10	10,000	5/VAC NAL-NL1 DSL v5.0a	✓	✓	Multiband Adaptive		✓			✓	✓	✓	✓	✓	✓	✓	*	✓	✓	✓		✓	MINI RITE / RITE BTE / MIC / CIC
Oticon	Agil	4	✓	✓	10	10,000	5/VAC NAL-NL1 DSL v5.0a	✓	✓	Multiband Adaptive		✓						✓			✓	*	✓	✓	✓		✓	MINI RITE / RITE BTE / MIC / CIC
Oticon	Dual XW	4 (Streamer)	✓	✓	10	10,000	Clarity 2 & 5/VAC	✓		Multiband Adaptive		✓			✓			✓	✓	✓	✓	312	Phonocoil				✓	
Oticon	Dual W	4 (Streamer)	✓	✓	8	8,000	Clarity 2 & 5/VAC	✓		Multiband Adaptive		✓								✓	✓	312	Phonocoil				✓	
Oticon	Dual V	4 (Streamer)	✓	✓	6	8,000	Clarity 2 & 5/VAC	✓		Adaptive		✓									✓	312	Phonocoil				✓	
Oticon	Vigo Pro connect	4		✓	8	8,000	3/NAL-NL1 & DSL 5.0	✓		Multiband Adaptive		✓									✓		VC ONLY				✓	RITE / BTE MIC / CIC
Oticon	Vigo connect	4		✓	6	8,000	3/NAL-NL1 & DSL 5.0	✓		Adaptive		✓									✓						✓	RITE / BTE MIC / CIC
Oticon	Dual m9			✓	10	10,000	Clarity 2 & 5/VAC	✓		Multiband Adaptive		✓						✓				10	Phonocoil				✓	
Oticon	Dual m7			✓	8	8,000	Clarity 2 & 5/VAC	✓		Multiband Adaptive		✓										10	Phonocoil				✓	
Oticon	Dual m5 / Pro			✓	6	8,000	Clarity 2 & 3/VAC	✓		Adaptive		✓										10/312	Phonocoil				✓	
Oticon	Hit Pro	4		✓	6	8,000	1 / NAL/NL1 & DSL v5.0a	✓		Adaptive		✓										*	✓	✓	✓		✓	BTE MIC / CIC
Oticon	Hit	4		✓	4	8,000	1 / NAL/NL1 & DSL v5.0a			Fixed		✓										*	✓	✓	✓		✓	BTE MIC / CIC
Oticon	Delta 4000				5	6,000	Clarity & 1/VAC			Fixed					DFC							10						

DIRECTIONALITY MODES	NOISE MANAGEMENT	FEEDBACK MANAGEMENT	BATTERY SIZE			
Multiband Adaptive Directionality Binaural Multiband Adaptive Directionality	Modulation Based TriState w/Voice Finder Binaural TriState w/Voice Finder	Dynamic Feedback Cancellation 2 Binaural Dynamic Feedback Cancellation 2	Battery size dependent upon model ordered	BINAURAL PROCESSING: Spatial Sound Binaural DFC My Voice	BINAURAL SYNCHRONISATION: Binaural Directionality Binaural Noise Management	BINAURAL COORDINATION: Binaural VC Binaural push button operation Binaural phone operation