

Product Information



Exélia Art

- 20 Channels
- Soundflow Premium
- SoundRecover
- 33 Band VoiceZoom
- NoiseBlock Processing
- WhistleBlock Technology
- WindBlock Manager
- EchoBlock System
- Datalogging+
- SoundRelax
- Selflearning
- Audiogram Direct
- NEW Insight
- Real Ear Sound
- DuoPhone with Direct Touch
- myPilot/iView*
- QuickSync
- ZoomControl with Direct Touch
- KeyPilot/WatchPilot
- iCom*
- iCube Programmable
- Easy Bluetooth
- Easy Audio
- Easy FM
- Easy Phone
- T-Coil
- AOV for custom products

* Dependent on model



Versáta Art

NEW

- 16 Channels
- Soundflow Advanced
- SoundRecover
- VoiceZoom
- NoiseBlock Processing
- WhistleBlock Technology
- WindBlock Manager
- Datalogging+
- User Preference Learning
- Audiogram Direct
- NEW Insight
- Real Ear Sound
- myPilot/iView*
- QuickSync
- KeyPilot2 / WatchPilot2
- Easy Bluetooth
- Easy Audio
- Easy FM
- Easy Phone
- T-Coil
- iCube Programmable
- AOV for custom products

* Dependent on model



Certéna Art

NEW

- 6 Channels
- Soundflow Standard
- SoundRecover
- Digital AudioZoom
- NoiseBlock Processing
- WhistleBlock Technology
- Datalogging+
- Audiogram Direct
- NEW Insight
- myPilot/iView
- QuickSync
- iCom
- KeyPilot2 / WatchPilot2
- Easy Bluetooth
- Easy Audio
- Easy Phone
- T-Coil
- iCube Programmable
- AOV for custom products



Milo Plus

- 4 Channels
- Digital AudioZoom
- NoiseBlock Processing
- WhistleBlock Technology
- Datalogging+
- Audiogram Direct
- T-Coil
- 4 manual accessible programs
- AOV available for custom products



Audéo SMART

- Up to 20 Channels*
- Soundflow Standard up to Premium*
- Digital AudioZoom up to VoiceZoom*
- SoundRecover
- SoundRelax†
- WindBlock Manager*
- EchoBlock System†
- DuoPhone† with Direct Touch
- Zoom Control† with Direct Touch
- Audiogram Direct
- New Insight
- Real Ear Sound*
- NoiseBlock Processing
- WhistleBlock Technology
- Datalogging+
- myPilot/iView
- QuickSync
- Push Button (for VC or Prog selection)
- Self Learning (IX only)
- iCom
- Easy Bluetooth
- Easy Audio
- Size 312 Battery
- Standard or Power Receiver
- iCube Programmable
- AOV on Custom Tips

* Dependent on model
† Audeo SMART IX only



Audéo MINI

- Up to 20 Channels*
- Soundflow Standard up to Premium*
- Digital AudioZoom up to VoiceZoom*
- SoundRecover
- SoundRelax†
- WindBlock Manager*
- EchoBlock System†
- Audiogram Direct
- New Insight
- NoiseBlock Processing
- WhistleBlock Technology
- Real Ear Sound
- Datalogging+
- Size 10 Battery
- Standard or Power Receiver
- AOV on Custom Tips

* Dependent on model
† Audeo MINI IX only

Product Information



Naída IX UP / SP

20 channels
4 manual accessible programs
Naida IX UP - 85/144 Ear sim (675)
Naida IX SP - 80/141 Ear sim (13)
VoiceZoom
Real Ear Sound
SoundRecover
ZoomControl with Direct Touch
Soundflow Premium
Bassboost
Easy FM
Quick Sync
WaterResistant
NoiseBlock
WhistleBlock
Audiogram Direct
DuoPhone with Direct Touch
EasyPhone
Easy Audio
Easy Bluetooth
iCom
iView
myPilot
T-Coil



Naída III UP / SP

6 channels
2+2 dedicated (FM and/or Tcoil)
Naida III UP - 85/144 Ear sim (675)
Naida III SP - 80/144 Ear sim (13)
AudioZoom
SoundRecover
Soundflow Standard
Bassboost
Easy FM
WaterResistant
NoiseBlock
WhistleBlock
Audiogram Direct
EasyPhone
Easy Audio
Easy Bluetooth
iCom
iView
myPilot
T-Coil

Features explained

Soundflow: A revolutionary system to give you the best possible hearing, wherever you are.

SoundRecover: High frequencies can be remapped to audible levels.

NoiseBlock Processing: Designed to maximise your comfort in noisy situations without affecting speech.

WhistleBlock Technology: The most effective anti-whistle system available today.

WindBlock Manager: A processor which will reduce the annoyance of wind turbulence.

EchoBlock: A system designed to help clarity wherever there are hard surfaces or high ceilings (eg Churches, rooms with laminate flooring).

myPilot: Remote control, alarm and program indicator, all with a full colour display.

iCom: Connect to bluetooth devices - mobile telephones, TV, SatNav.

Click'nTalk: Direct connection from selected mobile phones to any CORE wireless products without the need for iCom.

DuoPhone: Automatic bilateral hearing with reduced background noise with landline and mobile phones.

AudioGram Direct: A hearing test run through the instrument itself, to ensure your fitting is as personal as possible.

ZoomControl with Direct Touch: ZoomControl can now be set with the push of a button or using myPilot.



Exélia Art family with myPilot, iCom & Click'nTalk